

Solutran® hemo O₂O₃ Ozone therapy kit

by Meditalia

your solution

Kit for autohemotherapy with ozone

Graduated scale

To evaluate the quantity of the blood withdrawn.

Solutran

PVC phthalate free.

Swabable bag valve

Needle free.
To allow the connection of the transfusion set.

1-way valve

For the ozone addition without solution return.

Idrofobic filter

To allow the ozone addition.

Special connector for the ozone addition without using needle.

Technical notes:

- 450 ml phthalate free PVC bag
- Kit contains hermetically sealed automatic connectors
- Sterile with closed circuit
- **CE** 0373



Transfusion set

Flexible and kink resistant tubes.
Special components to guarantee safety, sterility and reliability of the kit.

Y connector

Accessible way for the addition of anticoagulant in case of blood clots in the part of the line between the needle and the Y connector.

Sodium citrate 3.8% vials



Solutran[®] hemo O₂O₃ Ozone therapy kit

by Meditalia

your solution

Properties

Sterile

The **Solutran O₂O₃ Kit** is “closed circuit”; the doctor can carry out the procedure in safety without the risk of environmental contamination.

Safe

The **Solutran O₂O₃ Kit** contains no needles, except the one used for drawing blood, throughout the circuit. The O₂O₃ mixture is added using a tube with no needle using special connectors.

No phthalate

The **Solutran O₂O₃ Kit** is phthalate free and is authorised for autohemotherapy with ozone.

Automatic connectors

The **Solutran O₂O₃ Kit** uses special antibacterial sealed connectors.

CE 0373 marked



Via Della Piana 1
23030 Lovero V. (SO) - Italy
meditaliasrl.com

Tel. +39 0342 77 10 70
Fax. +39 0342 77 10 71
info@meditaliasrl.com