

Solutran[®] hemo O₂O₃ Ozone therapy kit

by Meditalia

your solution

Properties

Sterile

The **Solutran O₂O₃ Kit** is “closed circuit”; the doctor can carry out the procedure in safety without the risk of environmental contamination.

Safe

The **Solutran O₂O₃ Kit** contains no needles, except the one used for drawing blood, throughout the circuit. The O₂O₃ mixture is added using a tube with no needle using special connectors.

No phthalate

The **Solutran O₂O₃ Kit** is phthalate free and is authorised for autohemotherapy with ozone.

Automatic connectors

The **Solutran O₂O₃ Kit** uses special antibacterial sealed connectors.

CE 0373 marked



VIA DELLA PIANA 1
23030 LOVERO V. (SO) - ITALY
MEDITALIASRL.COM

TEL. +39 0342 77 10 70
FAX. +39 0342 77 10 71
INFO@MEDITALIASRL.COM