


MEDICAL GRANULE



TRADE MARK Product			Solutran®	Solutran® Plus	Solutran® DOP Free	Solutran® Plus DOP Free	
Description			Medical Grade PVC Compound /with DEHP	Medical Grade CLEAR-PVC* Compound /with DEHP	Medical Grade PVC Compound - DEHP Free /with DEHT	Medical Grade CLEAR-PVC* Compound DEHP Free /with DEHT	
Application	Tubular Extrusion		●	●	●	●	
	Tube Extrusion		●		●		
Medical Use			For containers as a primary packaging for Infusion, Irrigation & Renal Systems	For containers as a primary packaging for Infusion, Irrigation & Renal Systems. For products that require exceptional transparency.	For containers as a primary packaging for Infusion, Irrigation & Renal Systems	For containers as a primary packaging for Infusion, Irrigation & Renal Systems. For products that require exceptional transparency.	
Plasticising	DEHP (DOP**)		●	●			
	DEHT (DOTP**)				●	●	
	TEHTM (TOTM**)						
	EVA						
Color ***	Natural		●		●		
	Natural Clear			●		●	
	Light Blue		●		●		
	Amber						
PHYSICAL PROPERTIES	Hardness	Value	60 Shore A up to 90 Shore A	80 Shore A up to 90 Shore A	60 Shore A up to 90 Shore A	80 Shore A up to 90 Shore A	
		Method	ISO 868	ISO 868	ISO 868	ISO 868	
	Melt Flow Index	Value ^(average)	-----	-----	-----	-----	
		Method	-----	-----	-----	-----	
	Density	Value ^(average)	1,2 g/cm ³	1,2 g/cm ³	1,2 g/cm ³	1,2 g/cm ³	
		Method	ISO R 1183	ISO R 1183	ISO R 1183	ISO R 1183	
	Tensile Breaking Load	Value ^(average)	From 135 to 170 Kg/cm ²	From 135 to 170 Kg/cm ²	From 130 to 175 Kg/cm ²	From 130 to 175 Kg/cm ²	
		Method	ISO R 527	ISO R 527	ISO R 527	ISO R 527	
	Elongation at Break	Value ^(average)	From 340% to 270% (depending on Shore)	From 340% to 270% (depending on Shore)	From 350% to 280% (depending on Shore)	From 340% to 270% (depending on Shore)	
		Method	ISO R 527	ISO R 527	ISO R 527	ISO R 527	
	Stiffening Temperature	Value ^(average)	From -22°C to -8°C (depending on Shore)	From -22°C to -8°C (depending on Shore)	From -22°C to -7°C (depending on Shore)	From -22°C to -8°C (depending on Shore)	
		Method	ISO R 458	ISO R 458	ISO R 458	ISO R 458	
	HANDLING	Storage Conditions	Packaging Type	Big Bags (Octagonal Containers) of 800/1000 Kg			
			Room Temperature	Not Exceeding 40 °C			
Temperature at Use			48h before use it should be stored within 18 °C to 22 °C				
Shelf-life		5 years from the date of production					

* Decrease of blushing effect after sterilization process

** Acronym In Italian Language

*** Color differences depend of the tone of the resin