Kifarm

your solution

by Meditalia



Technical notes:

- Capacity: 1000 cc
- Material: PE

400 00

300 cc

200 cc

100 00

- ETO sterilized
 - **C €** 0373 marked

Applications

Kifarm is to be used only by delivery room specializes staff member.

The intended use of Kifarm is to collect the bood loss during the childbirth and to allow a quickly visual estimation of it.







Kifarm

your solution

by Meditalia

Technical characteristics

• Thickness: 0.07 mm

• Lenght: 104 cm

• Width: 27 cm

 Scale with 100 cc graduation (from 0 to 1000 cc)

Packaging

- Primary packaging: medical paper printed blister
- Secondary packaging: box containing 200 bags individually packaged

